


IRANSCRIPTION

mRNA
PRODUCT
CATALOG
2024
V 2.4

Content

Schematic diagram of mRNA synthesis	01
Custom mRNA	02
BspQI (BSA-free)	03
High Yield T7 RNA Synthesis Kit	04
T7 RNA Polymerase	05
SP6 RNA Polymerase	06
Pyrophosphatase Inorganic (Yeast)	07
DNase I	08
RNase Inhibitor	08
Vaccinia Capping Enzyme	09
mRNA Cap 2'-O-Methyltransferase	10
Poly(A) polymerase	11
Other enzymes related to IVT	12
Nucleotides	13
Modified Nucleotides	13
Cap Analogs	13

Schematic diagram of mRNA synthesis



Custom mRNA

Good mRNA is far more than pure

Good mRNA should boost the translation after it is transfected or injected into cells, which requires a high integrity of the mRNA. In eukaryote, the cap1 structure at 5' end is vital because it promotes the translation, prevents the degradation and reduces the immunogenicity. The translation may be inhibited by the impurities and byproducts, such as the endotoxin, double stranded RNA, which are often be neglected.

Hongene mRNAs are produced under the control of GMP standard. Each mRNA is under a strict detection of the amount of dsRNA, capping efficiency, purity and integrity. So the mRNA is ready to be translated.


Feature

- Integrity > 90%
- Cap1 added to the 5' end, with the capping efficiency up to 99.8%
- 100 nt poly(A) added to the 3' end
- A260/A280 = 1.8-2.0
- A260/A230 > 2.0
- Endotoxin < 1 EU/mg RNA</p>
- dsRNA low to 0.002%
- Multiple choices for the modified nucleotides

Contact <u>info@hongene.com</u> to get the "Product Requirement Document for mRNA Synthesis Service", and your own mRNAs.

Alternative QC

mRNA sequencing mRNA size **Identity of Cap** Poly(A) tail length **Identity of modified bases Changes in base modification Identity of mRNA integrity** Residual NTP, SAM, SAH **Residual Enzymes Endotoxin Residual solvent** Bioburden Residual host cell protein **Residual host cell DNA Residual template DNA** Concentration **Capping efficiency** More...



Expression of Hongene eGFP.

2.5 µg eGFP mRNA was transfected into 1×10° HEK293T cell. 24 hours later, photos were taken under the fluorescent microscope.

Plasmid linearization — Ils restriction enzyme

BspQI (BSA-free) - No BSA

BspQI, which is cloned from *Bacillus sphaericus*, is a IIs restriction enzyme. It recognizes and cleaves at 5'-GCTCTTCN, $/N_c$ -3'.

5'- GCTCTTCNNNNN -3' 3'- CGAGAAGNNNNN -5'


It is recommended to linearize the template plasmid DNA at the site you like by a restriction enzyme to make sure the transcripts can be terminated accurately. The IIs restriction enzymes, such as BspQI is the most popular one because it does not introduce additional sequence at the cleavage site therefore gives a clean polyA tail. Meanwhile, the 3' overhang is avoided on the linearized template, which can reduce the amount of unexpected transcripts and the double stranded RNA byproducts.



Feature

 Animal-free, protein-free and endotoxin-free in reaction buffer.

Standability



Product information

Unit definition

One unit is defined as the amount of enzyme required to digest 1 μ g of λ DNA in 1X Cut Buffer following 1 hour incubation at 50°C in a total reaction volume of 50 μ L.

Concentration

10,000 units/mL

Reaction Temperature

50°C

Heat inactivation

80°C for 20 min

Storage buffer

20 mM Tris-HCl 1 mM DTT 0.1 mM EDTA 500 mM KCl 0.1% Triton X-100 50% Glycerol pH 7.0 @ 25°C

Reaction buffer

10X Cut Buffer (BSA-Free) 500 mM Tris-HCl 1000 mM NaCl 100 mM MgCl₂ pH 7.9 @ 25°C

Order information

Cat No.	Description	Amount	Storage
ON-124	BspQI	10 KU	-70°C/-20°C
N/A	10X Cut Buffer (BSA-Free)	5 mL	-20°C

Standability of BspQI of different concentrations under different storage conditions. The enzyme activities were assayed at 3, 5, 15, 19, 30, 34, 75, 151, 241, 357 days after the first storage.


RNA Synthesis—Kits High Yield T7 RNA Synthesis Kit

High Yield T7 RNA Synthesis Kit is designed to produce 25-fold more full-length RNA transcript per reaction than conventional in vitro transcription reactions. The yield is up to 210 µg RNA per reaction.

Usage

- Synthesize transcripts in vitro
- Synthesize gRNA for gene editing
- Synthesize a mRNA with fluorescent or radio-labeled nucleotides
- Synthesize a mRNA with base-modified nucleotides
- Synthesize a capped mRNA with cap analog in 1 step

Classic result




RNA synthesized with Hongene High Yield T7 RNA Synthesis Kit. The transcripts was cleanup by LiCl precipitation, and diluted by RNase-free water. 400 ng of the dilution was assayed with Agilent 5200 Fragment Analyzer.

Order information

Cat No.	Description	Amount	Storage
ON-040	High Yield T7 RNA Synthesis Kit	50 rxns	-20°C

Relevant Products

Cat No.	Description	Amount	Storage
R1331	ATP, 100 mM Sodium Solution	1 mL	-20°C
R2331	CTP, 100 mM Sodium Solution	1 mL	-20°C
R3331	GTP, 100 mM Sodium Solution	1 mL	-20°C
R5331	UTP, 100 mM Sodium Solution	1 mL	-20°C



RNA Synthesis—Enzyme T7 RNA Polymerase

T7 RNA polymerase is a monomeric bacteriophage encoded DNA directed RNA polymerase which catalyzes the formation of RNA in the 5' → 3' direction. In the process of initiation of transcription T7 recognizes a specific promoter sequence, the T7 promoter. T7 consists of 883 amino acids and has a molecular weight of 99 kDa.

T7 promoter

TAATACGACTCACTATAGGGAGA

Usage

- Synthesize the radio-labeled RNA probes
- Synthesize the RNA template for in vitro translation
- Research on RNA structure, process and catalytic mechanism

Feature

- Produce high yields in natural NTPs as well as modified ones
- Capping efficiency up to 95% when transcribing with some cap analogs
- Provide customized concentration to meet special requirements

Order information

Cat No.	Description	Amount	Storage
ON-004	T7 RNA polymerase	10 KU	-20°C
ON-062	10X IVT Reaction Buffer	5 mL	-20°C

Relevant Products

Cat No.	Description	Amount	Storage
R1331	ATP, 100 mM Sodium Solution	1 mL	-20°C
R2331	CTP, 100 mM Sodium Solution	1 mL	-20°C
R3331	GTP, 100 mM Sodium Solution	1 mL	-20°C
R5331	UTP, 100 mM Sodium Solution	1 mL	-20°C

Product information

Unit definition

One unit is defined as the amount of enzyme required to catalyze the incorporation of 1 nmol of ATP into acid-insoluble product in 1 hour at 37°C in a total volume of 20 µL.

Concentration

50,000 units/mL 200,000 units/mL 1000.000 units/mL


Storage buffer

20 mM Tris-HCl 100 mM NaCl 10 mM DTT 0.1% Triton X-100 1 mM EDTA 50%(v/v) Glycerol pH 7.9 @ 25°C

Reaction buffer

10X IVT Reaction Buffer 400 mM Tris-HCI 100 mM DTT 20 mM Spermidine 60 mM MgCl₂ pH 7.9 @ 25°C

Add 5'-cap by cotranscription



The transcriptional reaction without cap analog (upper) or with cap analog (down). After the reaction, the transcripts were annealed with 5' end DNA probes, digested by RNaseH. The purified 5' ends were assayed by HPLC. The capping efficiency was up to 95%.

Alternative choice for in vitro transcription —SP6 RNA Polymerase

SP6 RNA Polymerase is cloned from *Salmonella typhimurium* LT2Z, and expressed recombinantly. SP6 RNA Polymerase has a molecular weight of 98.5 kDa. It specifically recognizes sp6 promotor (as below), which is an alternative choice for gene design.

SP6 promoter +1 ATTTAGGTGACACTATAGAA

High Yield SP6 RNA Synthesis Kit is a kit used for in vitro transcription with SP6 RNA Polymerase. It is designed to produce 8-fold more full-length RNA transcript per reaction up to 80 µg with a reducing enzyme and NTPs consumption. Reagents for transcription, template depletion, purification and electrophoresis are concluded in the kit.

Usage

- Synthesize the radio-labeled RNA probes
- Synthesize the RNA template for in vitro translation
- Research on RNA structure, process and catalytic mechanism

Feature

- Produce high yields in natural NTPs as well as modified ones.
- Capping efficiency up to 95% when transcribing with some cap analogs.

Product information

Unit definition

One unit is defined as the amount of enzyme required to catalyze the incorporation of 1 nmol of ATP into acid-insoluble product in 1 hour at 37° C in a total volume of 20 µL.

Concentration

20,000 units/mL

Storage buffer

20 mM Tris-HCl 100 mM NaCl 20 mM DTT 0.1% Triton X-100 1 mM EDTA 50%(v/v) Glycerol pH 7.9 @ 25°C

Order information

Cat No.	Description	Amount	Storage
ON-338	High Yield SP6 RNA Synthesis Kit	50 rxns	-20°C
ON-297	SP6 RNA polymerase	2 KU	-20°C/-70°C
N/A	5X SP6 Reaction Buffer	0.5 mL	-20°C

RNA Synthesis—Enzyme


Pyrophosphatase Inorganic (Yeast)

Inorganic pyrophosphatase (Yeast) is prepared from an E. coli strain containing *Saccharomyces cerevisiae ppa* gene. It is a homodimer consisting of two equal subunits of molecular weight 32-35 kDa. Inorganic pyrophosphatase catalyzes the following reaction:

$P_2O_7^{-4} + H_2O \rightarrow 2HPO_4^{-2}$

The synthesis of RNA will generate pyrophosphate, which can bind Mg²⁺ to form sediment and result in a declined yield of RNA. However, under the action of pyrophosphatase, the pyrophosphate byproduct can be degraded into phosphate, which pulls the reaction far in the synthetic direction.

IPPA influence the yield of transcription



Inorganic pyrophosphatase (Yeast) was added into a 20 μ L IVT reaction at a final concentration of 0.0125-1 U/mL. The mRNA yield was enhanced.

Product information

Unit definition

One unit is the amount of enzyme that will generate 1 µmol of phosphate per minute from inorganic pyrophosphate under standard reaction conditions.

Concentration

100 units/mL

Storage buffer

20 mM Tris-HCl 1 mM DTT 0.1 mM EDTA 50% Glycerol pH 8.0 @ 25°C

Order information

Cat No.	Description	Amount	Storage
ON-025	Pyrophosphatase Inorganic (Yeast)	10 KU	-20°C

Relevant Products

Cat No.	Description	Amount	Storage
ON-004	T7 RNA polymerase	10 KU	-20°C
ON-062	10X IVT Reaction Buffer	0.5 mL	-20°C
ON-039	Ribonuclease Inhibitor, Human Placenta	40 KU	-20°C

Template Digestion --DNase I

DNase I is an endonuclease that digests singleand double-stranded DNA. It nonspecifically cleaves DNA to release di-,tri- and oligonucleotide products with 5'-phosphorylated and 3'-hydroxylated end.

Usage

Digest the template DNA after transcriptional reaction

Order information

Cat No.	Description	Amount	Storage
ON-109	DNase I (RNase-free)	1000 U	-20°C
ON-077	10X DNase I Reaction Buffer	1 mL	-20°C

Product information

Unit definition

One unit is defined as the amount of enzyme required to completely degrade 1 μ g of λ DNA in 10 minutes at 37°C in a total volume of 50 μ L.

Concentration

1,000 units/mL

Storage buffer 50 mM Tris-HCl 10 mM CaCl₂ 50% Glycerol

pH 7.5 @ 25°C

Reaction buffer

10X DNase I Reaction Buffer 100 mM Tris-HCI 25 mM MgCl₂ 10 mM CaCl₂ pH 7.6 @ 25°C

Prevent Degradation —RNase Inhibitor

Product information

Unit definition

One unit is defined as the amount of Recombinant Ribonuclease Inhibitor required to inhibit the activity of 5 ng of ribonuclease A by 50%. Activity is measured by the inhibition of hydrolysis of cytidine 2'3'-cyclic monophosphate by ribonuclease A.

Concentration

40,000 units/mL

Storage buffer

20 mM HEPES-NaOH 50 mM NaCl 8 mM DTT 50%(v/v) glycerol pH 7.6 @ 25°C **Ribonuclease Inhibitor, Human Placenta** is a recombinant human placental protein with specifically inhibits broad-spectrum RNase such as RNase A, RNase B and RNase C. But It is not effective against RNase I, RNase TI, SI Nuclease, RNase H, Taq DNA polymerase, M-MLV Reverse Transcriptase and T7 RNA Polymerase. The 50kDa protein exerts its inhibitory effect by noncovalently binding to RNases at 1:1 ratio. The Ki value for binding of RNase Inhibitor, Human Placenta to RNase is approximately 10-14 M.

<u>Usage</u>

- In vitro transcription/translation
- cDNA synthesis
- RNA labeling by enzymatic method

Order information

Cat No.	Description	Amount	Storage
ON-039	Ribonuclease Inhibitor, Human Placenta	40 KU	-20°C


Capping —Vaccinia Capping Enzyme

Vaccinia Capping Enzyme is a recombinant protein from Vaccinia virus. Vaccinia Capping Enzyme can add 7-methylguanylate cap structure (Cap 0) to the 5' end of RNA.

In eukaryotes, these terminal cap structures are involved in stabilization, transport and translation of mRNAs.

RNA triphosphatase, RNA quanylyltransferase, and RNA (guanine-7)-methyltransferase activities are associated with the vaccinia virus RNA capping enzyme, a heterodimeric protein containing polypeptides of Mr 95,000 and Mr 31,000.

In vitro, transcripts can be capped in the presence of reaction buffer, GTP, and the methyl donor (SAM).



Add 5'-cap by Vaccinia Capping Enzyme

Usage

- Add natural cap structure to 5' end of mRNA
- Label the 5' end of mRNA
- Synthesize Cap1 mRNA together with 2'-O-Methyltransferase

Product information

Unit definition

One unit is defined as the amount of enzyme required to incorporate 10pmol of GTP into an 80 nt transcript in 1 hour at 37°C.

Concentration

10,000 units/mL

Storage buffer 20 mM Tris-HCl 0.1 mM EDTA 100 mM NaCl 1 mM DTT 0.1%(v/v) Triton X-100 10 mM DTT 50% Glycerol pH 8.0 @ 25°C

Reaction buffer

10X Capping Buffer 500 mM Tris-HCl 50 mM KCl 10 mM MgCl₂ pH 8.0 @ 25°C

Brief protocol

- 1. Heat at 65°C for 5 minutes, Ice bath for 5 minutes.
- 2. Prepare a 20 µL reaction, 2'-O-Methyltransferase can be added simultaneously.
- 3. Incubation at 37°C for 30 minutes.

Order information

Cat No.	Description	Amount	Storage
ON-028	Vaccinia Capping Enzyme	10 KU	-20°C
ON-073	10X Capping Buffer	2 mL	-20°C
ON-074	S-adenosylmethionine (SAM)	1 mL	-20°C
ON-075	10 mM GTP	1 mL	-20°C

Relevant Products

Cat No.	Description	Amount	Storage
ON-014	mRNA Cap 2'-O-Methyl Transferase	50 KU	-20°C

Capping

--mRNA Cap 2'-O-Methyl-transferase

mRNA Cap-2'-O-Methyltransferase is a recombinant protein from Vaccinia virus. The enzyme can add a methyl group at the 2'-O position of the first nucleotide adjacent to the Cap structure at the 5' end of the RNA. The enzyme utilizes SAM as a methyl donor to methylate capped RNA (Cap 0) resulting in a Cap 1 structure.

The Cap structure has been reported to enhance mRNA translation efficiency and hence may help improve expression in mRNA transfection and microinjection experiments.

mRNA cap-2'-O-Methyltransferase specifically requires RNA with an m7GpppN Cap as substrate.

Usage

- Transfer Cap0 mRNA into Cap1 mRNA
- Enhance translation efficiency of mRNA

Order information

Cat No.	Description	Amount	Storage
ON-014	mRNA Cap 2'-O- Methyltransferase	50 KU	-20°C
ON-073	10X Capping Buffer	2 mL	-20°C
ON-074	S-adenosylmethionine (SAM)	1 mL	-20°C

Product information

Unit definition

One unit is defined as the amount of enzyme required to methylate 10 pmols of 80 nts Capped RNA transcript in 1 hour at 37°C.

Concentration

50,000 units/mL

Storage buffer

20 mM Tris-HCl
0.1 mM EDTA
100 mM NaCl
1 mM DTT
0.1%(v/v) Triton X-100
50% Glycerol
pH 8.0 @ 25°C

Reaction buffer

10X Capping Buffer 500 mM Tris-HCl 50 mM KCl 10 mM MgCl₂ 10 mM DTT pH 8.0 @ 25°C

Tailing

--Poly(A) polymerase

Poly(A) Polymerase uses ATP as a substrate for template-independent addition of adenosine monophosphate to the 3'-hydroxyl termini of RNA molecules.

Usage

- Add polyA tail to transcripts
- Enhance the translation efficiency
- Label RNA with fluorescent or radio-labeled ATP

Order information

Cat No.	Description	Amount	Storage
ON-126	Poly(A) Polymerase	5000 U	-20°C
ON-127	10X poly(A) Polymerase Reaction Buffer	2 mL	-20°C
R1331	ATP, 100mM Sodium Solution	1 mL	-20°C

Relevant products

Cat No.	Description	Amount	Storage
ON-040	High Yield T7 RNA Synthesis Kit	50 rxns	-20°C
ON-039	Ribonuclease Inhibitor, Human Placenta	40 KU	-20°C

Product information

Unit definition

One unit is defined as the amount of enzyme that will incorporate 1 nmol of AMP into RNA in a 20 μ L volume in 10 min at 37°C.

Concentration

5,000 units/mL

Storage buffer

20 mM Tris-HCl 1 mM EDTA 300 mM NaCl 1 mM DTT 0.1%(v/v) Triton X-100 50% Glycerol pH 7.5 @ 25°C

Reaction buffer

10X poly(A) Polymerase Reaction Buffer 500 mM Tris-HCl 2.5 M NaCl 100 mM MgCl₂ pH 7.9 @ 25°C

RNA Synthesis

Nucleotides and Cap Analogs

Hongene Biotech Corporation has become a manufacturer with state-of-art facilities in the fields of nucleosides, nucleotides and phosphoramidites since 2001.

At Hongene our experts are also happy to offer custom synthesis products essential for your R&D, helping to bring your products to market more efficiently. Through continuous innovation, improvement of manufacturing technologies, the tightest quality control and devoted technical team and services, we aim to help our customers to achieve their goals efficiently and economically.

NTPs

To customize cap analogs, please contact info@hongene.com

Cat No.	Description	Amount	Storage
R1331	ATP, 100 mM Sodium Solution	1 mL	-20°C
R2331	GTP, 100 mM Sodium Solution	1 mL	-20°C
R3331	CTP, 100 mM Sodium Solution	1 mL	-20°C
R5331	UTP, 100 mM Sodium Solution	1 mL	-20°C
R5-027	N1-Me- Pseudo UTP, 100 mM Sodium Solution	100 µL	-20°C
R1-051	ATP, 200 mM Tris Solution	100 µL	-20°C
R2-057	GTP, 200 mM Tris Solution	100 µL	-20°C
R3-052	CTP, 200 mM Tris Solution	100 µL	-20°C
R5-065	UTP, 200 mM Tris Solution	100 µL	-20°C
R5-064	N1-Me-Pseudo UTP, 200mM Tris Solution	100 µL	-20°C
R3-029	5-Me-CTP, 100mM Sodium Solution	1 mL	-20°C
R5-104	5-Me-UTP, 100mM Sodium Solution	1 mL	-20°C
R5-022	Pseudo UTP, 100mM Sodium Solution	1 mL	-20°C
R5-046	5-OMe-UTP, 100mM Sodium Solution	1 mL	-20°C
R5-066	5-OMe-UTP, 200mM Tris Solution	100 µL	-20°C

Cap Analogs

Cat No.	Description	Amount	Storage
ON-089	GpppG	Inquiry	-20°C
ON-137	GpppA	Inquiry	-20°C
ON-138	M7-GpppA	Inquiry	-20°C
ON-136	M7-GpppG	Inquiry	-20°C
ON-134	GAG Trimer*	Inquiry	-20°C

^{*}The use of this product may require the buyer to obtain additional third party intellectual property rights for certain applications.

Enzymes & Kits

Kits

Cat No.	Description	Amount	Storage
ON-040	High Yield T7 RNA Synthesis Kit	50 rxns	-20°C
ON-257	Capped RNA Synthesis Kit*	N/A	-20°C
ON-269	Capped RNA Synthesis Kit (with tailing)*	N/A	-20°C

^{*} For sale outside USA only.

Enzymes

Cat No.	Description	Amount	Storage
ON-124	BspQI (BSA-free)	10 KU	-70°C/-20°C
ON-211	T7 Enzyme mix	N/A	-20°C
ON-004	T7 RNA Polymerase, low concentration	10 KU	-20°C
ON-005	T7 RNA Polymerase, high concentration	N/A	-20°C
ON-039	RNase Inhibitor (recombinant)	40 KU	-20°C
ON-025	Pyrophosphatase, Inorganic (yeast)	10 KU	-20°C
ON-109	DNase I (recombinant, RNase-free)	1000 U	-20°C
ON-028	Vaccinia Capping Enzyme	10 KU	-20°C
ON-014	2'-O-Methyltransferase	50 KU	-20°C
ON-074	S-adenosylmethionine (SAM), 32mM	N/A	-20°C
ON-126	Poly(A) Polymerase	5000 U	-20°C
ON-024	RNase III, <i>E.coli</i>	1000 U	-70°C
ON-078	10X RNase III Reaction Buffer	1 mL	-20°C
ON-079	10X MnCl ₂	1 mL	-20°C
ON-080	10X EDTA	1 mL	-20°C
ON-081	T4 RNA Ligase 1	N/A	-20°C
ON-179	Alkaline Phosphatase, TAB5	1000 U	-20°C
ON-180	10X Phosphatase Reaction Buffer	0.8 mL	-20°C
ON-333	NLS-Cas9, 10 µM	1000 pmol	-70°C

CDMO Services Available



HEADQUARTER

Shanghai Hongene Biotech Corporation 9F, Building 1#, No.315 Guiping Road Shanghai, China, 200233

Tel: +86-21-6475-7213 Fax: +86-21-6470-0613 Email: sales@hongene.com

JAPAN OFFICE

Hongene Biotech Japan Co., Ltd. 403, Kanda Iwamotocho Plaza Bldg, 2-4-1 Iwamotocho,

Chiyoda-ku, Tokyo 101-0032, Japan

Tel: +81-3-5829-9770 Fax: +81-3-5829-9795 Email: info@hongene.jp

US OFFICE

Hongene Biotech Corporation 29520 Kohoutek Way Union city CA.94587

Tel: +1-510-931-4711

Primary contact: product.info@hongene.com Secondary contact: info@hongene.com

GERMANY OFFICE

Hongene Biotech Germany GmbH Großmoorbogen 25, 21079 Hamburg, Germany

Phone: +49 40 8830 6130

Primary contact: info-eu@hongene.com Secondary contact: info@hongene.com

> GMP Facility β-Lactams-Antibiotics-FREE Animal-FREE